


## CERTIFICATE OF ANALYSIS

<b>Product Name:</b> 2-Chloro-N-ethylethanamine hydrochloride		
<b>Chemical Name:</b> N-Ethyl-N-(2-chloroethyl)amine hydrochloride		
<b>CAS number</b>	4535-87-9	<b>Chemical structure</b>  
<b>Molecular Formula</b>	C4H11Cl2N	
<b>Molecular Weight</b>	144.04	
<b>Analysis Date</b>	March-2025	
<b>Re-test Date</b>	March-2027	
<b>Batch Number</b>	VLS-08-CEE-137	
<b>Storage condition</b>	Store at 2-8°C	
<b>Shipment condition</b>	At Room temperature	

### ANALYTICAL INFORMATION:

S No	Test	Results
1	Description	Off-white color hygroscopic solid (should be handled under N <sub>2</sub> )
2	Solubility	H <sub>2</sub> O, DMSO (Slightly soluble)
3	Identification <div>             1. <sup>1</sup>H-NMR              2. Mass           </div>	Confirms to the structure  Confirms to the structure
4	Chromatographic purity (By HPLC)	95.11%

Prepared by:	Reviewed by:
--------------	--------------