


## CERTIFICATE OF ANALYSIS

<b>Product Name: 2,3-dichloro-4-nitro-benzenemethanamine</b>		
<b>Chemical Name: 2,3-dichloro-para-nitro benzyl amine</b>		
<b>CAS number</b>	1806356-69-3	<b>Chemical structure</b>  
<b>Molecular Formula</b>	C <sub>7</sub> H <sub>6</sub> Cl <sub>2</sub> N <sub>2</sub> O <sub>2</sub>	
<b>Molecular Weight</b>	221.04	
<b>Analysis Date</b>	March-2025	
<b>Re-test Date</b>	March-2027	
<b>Batch Number</b>	VLS-08-ANG-121	
<b>Storage condition</b>	Store at 2-8°C	
<b>Shipment condition</b>	At Room temperature	

### ANALYTICAL INFORMATION:

S No	Test	Results
1	<b>Description</b>	Brown color Hygroscopic solid
2	<b>Solubility</b>	DMSO
3	<b>Identification</b>  1. 1H-NMR  2. Mass	Conforms to the structure  Conforms to the structure
4	<b>Chromatographic purity (By HPLC)</b>	<b>93.12%</b>

Prepared by:	Reviewed by:
--------------	--------------