


CERTIFICATE OF ANALYSIS

Product Name : (Z)-5-fluoro-3-hydrazonoindolin-2-one

Chemical Name: (Z)-5-Fluoro-1H-indole-2,3-dione-3-Hydrazone

CAS number	283584-52-1	Chemical structure 
Molecular Formula	C ₈ H ₆ FN ₃ O	
Molecular Weight	179.15	
Analysis Date	January-2025	
Re-test Date	January-2027	
Batch Number	VLS-06/SAA/05	
Storage condition	At 2-8°C	
Shipment condition	At room temperature	

ANALYTICAL INFORMATION:

S No	TEST	RESULT
1	DESCRIPTION	Yellow color solid
2	SOLUBILITY	DMSO
2	IDENTIFICATION	
	1.H ¹ -NMR 2.MASS	Conforms to the structure Conforms to the structure
3	CHROMATOGRAPHIC PURITY (By HPLC), area Normalization)	91.67%

Prepared by:

Reviewed by: