


CERTIFICATE OF ANALYSIS

Product Name : Abiraterone Isopropyl Ether

Chemical Name: 3beta-Isopropoxy-17-(pyridin-3-yl)androsta-5,16-diene

CAS number	2484719-15-3	Chemical structure 
Molecular Formula	C27H37NO	
Molecular Weight	391.60	
Analysis Date	January-2025	
Re-test Date	January-2027	
Batch Number	VLS-04/ARB-I/240	
Storage condition	At 2-8°C	
Shipment condition	At room temperature	

ANALYTICAL INFORMATION:

S No	TEST	RESULT
1	DESCRIPTION	White semi solid powder
2	SOLUBILITY	DMSO
2	IDENTIFICATION 1.H ¹ -NMR 2.MASS 3.IR	Conforms to the structure Conforms to the structure Conforms to the structure
3	CHROMATOGRAPHIC PURITY (By HPLC)	95.58%

Prepared by:

Reviewed by: