


## CERTIFICATE OF ANALYSIS

Product Name : 3-(DIMETHYLAMINO)-1-(NAPHTHALEN-2-YL)PROPAN-1-OL

Chemical Name: 3-(DIMETHYLAMINO)-1-(NAPHTHALEN-2-YL)PROPAN-1-OL

CAS number	13634-66-7	<b>Chemical structure</b>  
Molecular Formula	C <sub>15</sub> H <sub>19</sub> NO	
Molecular Weight	229.32	
Analysis Date	November-2024	
Re-test Date	November-2026	
Batch Number	VLS-04/BDQ/11	
Storage condition	At -20°C	
Shipment condition	At Room temperature	

### ANALYTICAL INFORMATION:

S No	TEST	RESULT
1	DESCRIPTION	White solid
2	Solubility	DMSO, MeOH
2	IDENTIFICATION 1.H <sup>1</sup> -NMR 2.Mass 3.IR	Conforms to the structure Conforms to the structure Conforms to the structure
3	CHROMATOGRAPHIC PURITY (By HPLC)	98.46%

Prepared by:

Reviewed by: