


CERTIFICATE OF ANALYSIS

Product Name : (E)-N-(4-fluorophenyl)-2-(hydroxyimino)acetamide

Chemical Name: (E)-N-(4-fluorophenyl)-2-(hydroxyimino)acetamide

CAS number	879274-72-3	Chemical structure 
Molecular Formula	C ₈ H ₇ FN ₂ O ₂	
Molecular Weight	182.15	
Analysis Date	January-2025	
Re-test Date	January-2027	
Batch Number	VLS-05/SAA/229	
Storage condition	At 2-8°C	
Shipment condition	At room temperature	

ANALYTICAL INFORMATION:

S No	TEST	RESULT
1	DESCRIPTION	Grey color solid
2	SOLUBILITY	DMSO
2	IDENTIFICATION 1.H ¹ -NMR 2.MASS	Conforms to the structure Conforms to the structure
3	CHROMATOGRAPHIC PURITY (By HPLC), area Normalization)	93.12%

Prepared by:

Reviewed by: