


CERTIFICATE OF ANALYSIS

Product Name: EZE-5 Nitroso

Chemical Name: N-(4-Fluorophenyl)-N-((1S,2R,5S)-5-(4-fluorophenyl)-5-hydroxy-1-(4-hydroxyphenyl)-2-((S)-2-oxo-4-phenyloxazolidine-3-carbonyl)pentyl)nitrous Amide

CAS number	NA	Chemical structure 
Molecular Formula	C ₃₃ H ₂₉ F ₂ N ₃ O ₆	
Molecular Weight	601.6	
Analysis Date	March-2025	
Re-test Date	March-2027	
Batch Number	VLS-08/EZE/157	
Storage condition	Store at 2-8°C	
Shipment condition	At Room temperature	

ANALYTICAL INFORMATION:

S No	Test	Results
1	Description	Yellow color solid
2	Solubility	DMSO
3	Identification	
	1. 1H-NMR	Conforms to the structure
	2. Mass	Conforms to the structure
4	Chromatographic purity (By HPLC)	95.11%

Prepared by:

Reviewed by: