


CERTIFICATE OF ANALYSIS

Product Name : *N*-Phenylacetamide

Chemical Name: Acetanilide

CAS number	103-84-4	Chemical structure 
Molecular Formula	CH ₃ CONHC ₆ H ₅	
Molecular Weight	135.16	
Analysis Date	January-2025	
Re-test Date	January-2027	
Batch Number	VLS-04/NPA/243	
Storage condition	At 2-8°C	
Shipment condition	At room temperature	

ANALYTICAL INFORMATION:

S No	TEST	RESULT
1	DESCRIPTION	White semi solid powder
2	SOLUBILITY	DMSO
2	IDENTIFICATION 1.H ¹ -NMR 2.MASS 3.IR	Conforms to the structure Conforms to the structure Conforms to the structure
3	CHROMATOGRAPHIC PURITY (By HPLC)	95.58%

Prepared by:	Reviewed by:
--------------	--------------