


## CERTIFICATE OF ANALYSIS

<b>Product Name:</b> EZE-5 Impurity A Nitroso		
<b>Chemical Name:</b> (S)-(1-nitrosopyrrolidin-2-yl)diphenylmethanol		
<b>CAS number</b>	NA	<b>Chemical structure</b>  
<b>Molecular Formula</b>	C <sub>17</sub> H <sub>18</sub> N <sub>2</sub> O <sub>2</sub>	
<b>Molecular Weight</b>	282.34	
<b>Analysis Date</b>	March-2025	
<b>Re-test Date</b>	March-2027	
<b>Batch Number</b>	VLS-08/EZE/153	
<b>Storage condition</b>	Store at 2-8°C	
<b>Shipment condition</b>	At Room temperature	

### ANALYTICAL INFORMATION:

S No	Test	Results
1	Description	Off-white color solid
2	Solubility	DMSO
3	Identification 1. 1H-NMR 2. Mass	Conforms to the structure Conforms to the structure
4	Chromatographic purity (By HPLC)	99.97%

Prepared by:	Reviewed by:
--------------	--------------