


CERTIFICATE OF ANALYSIS

Product Name : Tenofovir related compound B

Chemical Name: (E)-9-(Prop-1-enyl)-9H-purin-6-amine

CAS number	1446486-33-4	Chemical structure 
Molecular Formula	C ₈ H ₉ N ₅	
Molecular Weight	175.19	
Analysis Date	November-2024	
Re-test Date	November -2026	
Batch Number	VLS-04/TAA/97	
Storage condition	At -20°C	
Shipment condition	At room temperature	

ANALYTICAL INFORMATION:

S No	TEST	RESULT
1	DESCRIPTION	White color solid
2	SOLUBILITY	DMSO
2	IDENTIFICATION 1.H ¹ -NMR 2.MASS 3.IR	Conforms to the structure Conforms to the structure Conforms to the structure
3	CHROMATOGRAPHIC PURITY (By HPLC), area Normalization)	93.13%

Prepared by:

Reviewed by: