


CERTIFICATE OF ANALYSIS

Product Name : ERLOTINIB

Chemical Name: 2-METHOXY ETHYL 2-AMINO 4,5-BIS (2-METHOXY ETHOXY) BENZOATE

CAS number	NA	Chemical structure 
Molecular Formula	C ₁₆ H ₂₅ NO ₇	
Molecular Weight	343.38	
Analysis Date	December-2024	
Re-test Date	December-2026	
Batch Number	VLS-05/ERL/29	
Storage condition	At -20°C	
Shipment condition	At Room temperature	

ANALYTICAL INFORMATION:

S No	TEST	RESULT
1	DESCRIPTION	Brown color residue
2	SOLUBILITY	DMSO
2	IDENTIFICATION 1.H ¹ -NMR 2.MASS 3.IR	Conforms to the structure Conforms to the structure Conforms to the structure
3	CHROMATOGRAPHIC PURITY (By HPLC), area Normalization)	96.62%

Prepared by:

Reviewed by: