


CERTIFICATE OF ANALYSIS

Product Name : N-(Chloro(dimethylamino)methylene)-N-methylmethanaminium hexafluorophosphate(V)

Chemical Name : Chloro-N,N,N',N'-tetramethylformamidine Hexafluorophosphate

CAS number	94790-35-9	Chemical structure 
Molecular Formula	C ₅ H ₁₂ ClF ₆ N ₂ P	
Molecular Weight	280.58	
Analysis Date	January-2025	
Re-test Date	January-2027	
Batch Number	VLS-05/INT-2/169	
Storage condition	At -20°C	
Shipment condition	At room temperature	

ANALYTICAL INFORMATION:

S No	TEST	RESULT
1	DESCRIPTION	White Color solid
2	SOLUBILITY	DMSO
3	IDENTIFICATION 1.H ¹ -NMR 2.MASS 3.IR	Conforms to the structure Conforms to the structure Conforms to the structure
4	CHROMATOGRAPHIC PURITY (By HPLC), area Normalization)	99.97%

Prepared by:

Reviewed by: