


CERTIFICATE OF ANALYSIS

Product Name: N-Nitroso Amlodipine

Chemical Name: 3-ethyl 5-methyl 2-((2-aminoethoxy)methyl)-4-(2-chlorophenyl)-6-methyl-1-nitroso-1,4-dihydropyridine-3,5-dicarboxylate

CAS number	NA	Chemical structure 
Molecular Formula	C ₂₀ H ₂₄ ClN ₃ O ₆	
Molecular Weight	437.88	
Analysis Date	March-2025	
Re-test Date	March-2027	
Batch Number	VLS-08-ANG-121	
Storage condition	Store at 2-8°C	
Shipment condition	At Room temperature	

ANALYTICAL INFORMATION:

S No	Test	Results
1	Description	Yellow color solid
2	Solubility	DMSO
3	Identification 1. 1H-NMR 2. Mass	Conforms to the structure Conforms to the structure
4	Chromatographic purity (By HPLC)	95.12%

Prepared by:

Reviewed by: