


CERTIFICATE OF ANALYSIS

Product Name : 3-(5-Chloro-4-(pyridin-3-yl)pyrimidin-2-ylamino)-4-methyl-N-(3-(4-methyl-1H-imidazol-1-yl)-5-(trifluoromethyl)phenyl)benzamide

Chemical Name: 3-(5-Chloro-4-(pyridin-3-yl)pyrimidin-2-ylamino)-4-methyl-N-(3-(4-methyl-1H-imidazol-1-yl)-5-(trifluoromethyl)phenyl)benzamide

CAS number	2469039-94-7	Chemical structure 
Molecular Formula	C ₂₆ H ₂₁ ClF ₃ N ₇ O	
Molecular Weight	563.96	
Analysis Date	May-2025	
Re-test Date	April-2027	
Batch Number	VLS-09-NTB-CL-B-29	
Storage condition	At 2-8°C	
Shipment condition	At room temperature	

ANALYTICAL INFORMATION:

S No	Test	Results
1	Description	Pale pink to beige color solid
2	Solubility	DMSO(Slightly)
3	Identification	
	1. 1H-NMR	Conforms to the structure
	2. Mass	Conforms to the structure
4	Chromatographic purity (By HPLC)	91.18%

Prepared by:	Reviewed by:
--------------	--------------